UK

# Clippo bedwetting alarm system

Product information

User manual

Warranty

**rodger**

## **Thanks, you've selected the Rodger bedwetting alarm**

Thank you for your confidence in the Rodger bedwetting alarm. We assume that, if you use the bedwetting alarm in a good way, you will have success. Please, read all instructions before using the bedwetting alarm.

- The Rodger bedwetting system is suitable for children over the age of 5 years.
- When beginning the therapy, be motivated, relaxed and confident.
- Test the alarm before the start of the treatment. You and your child will get familiar with the bedwetting alarm. Repeat this before your child goes to sleep.
- Make no contest of the therapy. This will cause stress and is counterproductive.
- Give your child an occasional reward. A small present is a nice positive reinforcement that your child is doing well.
- Use the Rodger bedwetting system every night until there have been no "wet beds" for 14 consecutive nights.
- Successful treatment with the Rodger bedwetting system takes from a few weeks to a few months. Persistence is the key to success.
- Do your best to use the alarm every night without interruption.
- If your child's bedwetting condition does not improve within 4 months, visit your doctor to rule out any medical conditions.
- Do not start the therapy if important events will take place soon, i.e. moving the family or the birth of a child.
- Do not wake your child at night to go to the toilet.
- Never punish your child; your child is not responsible for occasional accidents.
- Do not prohibit your child from drinking before going to bed. Drinking less does not decrease the likelihood of bedwetting. However, it is preferable that your child does not drink coffee, tea, or beverages containing carbonation or chocolate before going to bed.
- In the early days of the treatment, a child who is a heavy sleeper might not wake up. You could choose another sound of the alarm.

Another good solution could be to sleep next to your child and help your child to wake up as soon as the alarm sounds, but make sure that your child stops the alarm and goes to the toilet. In time, your child will “tune in” and do this more independently.

- If relapse occurs, another short course will usually succeed. Store all parts together for future use. Remove battery before storing.

## Explanation Rodger Clippo

See inside cover for an overview of the content.


- A Sensor cable connection
- B Stop button
- C Alarm type switch
- D Magnetic clip
- E Sound type switch
- F Volume switch
- G Select button
- H Sensorclip
- I Rodger pantsconnector (optional)
- J Rodger pants (optional)

## Settings

See inside cover for detailed pictures.

- **Setting alarm type:**  
Select the alarm type of your choice with button (C):
  1. Just vibration function (V)
  2. Sound and vibration function (V&S)
  3. Just sound function (S)
- **Setting sound type (remove battery cover):**  
**Option 1:** For every alarm to be the same sound > set switch (E) to 1;  
Select the sound of your choice by pressing the “select” button (G).  
The last selected sound will appear every time when the alarm goes off.

To test the sound, put the “test device” in the mouth of the sensorclip. This will cause the alarm to sound. (In case you ordered the optional pantsconnector, connect the two push-buttons of the pantsconnector with the “test device”).

**Option 2:** For every alarm to be a different sound > set switch **(E)** to 8;

With every alarm 1 of the 8 available sounds will sound.

- **Setting volume (remove battery cover):**  
Select the desired volume with switch **(F)**.

## Manual instruction pictograms

See inside cover for an overview of the instruction pictograms belonging to the explanations below.

1. Connect the cable to the bedwetting alarm.
2. Clip the alarm with the magnetic clip to the shoulder area of the pajama top.
3. Attach the sensorclip to the outside of the pants.
4. Sleep well.
5. The alarm will sound when urination occurs.
6. Get up, remove the clip from the pants, turn off the alarm and go to the toilet.
7. Repeat steps 2 and 3 with a clean pants, sleep well.
8. Wake up dry!!!


## Remarks

- Attach the sensorclip in the area that the first drops of urine would appear (see figure).
- Let the cable run inside the pajama top.
- Only removing the clip will not stop the alarm.
- Clean the sensorclip regularly. You can do this by using a little soapy water on a cloth. Rinse afterwards with clean water and dry.

## Optional: Rodger pantsconnector and Rodger pants

For a more accurate detection you can order the optional Rodger pantsconnector in combination with special Rodger pants. These pants are equipped with sensor wires for immediate detection of urine. See inside cover for images of the Rodger pantsconnector **(I)** and the Rodger pants **(J)**.

Simply click the pantsconnector on the pants.

### Laundering Rodger pants:

The pants need to be washed before use. They become more absorbent. The pants can be machine washed and tumble dried. Do certainly not wring! Avoid fabric softener. Over time, residue may build up on the sensor contacts inside the pants. You simply need to soak the sensor pants for about 30 minutes in vinegar solution (1 cup vinegar and 2 cups water).

### Testing Rodger pants:

- Take a cup of lukewarm saltwater (1 teaspoon salt in 1 cup water).
- Click the pantsconnector on the pants.
- Connect the pantsconnector to the system.
- Pour two teaspoons of this salty solution on the front of the pants.  
This should cause the alarm to sound.
- Remove the pantsconnector from the pants.
- The alarm can be turned off by pushing the stop button **(B)**.

## Online success calendar

The online success calendar is an interactive version of the standard success calendar and is very motivating during the treatment. By using the online success calendar, you receive informative messages about the progress of the training. After a dry night, your child may choose a game of its choice as a reward. Further instructions: [www.rodger.nl](http://www.rodger.nl).


## Explanation success calendar

	1	2		4			
							
							
							
							

Dry: Slept through

Dry: Visited toilet

Wet: Woke by alarm

Wet: Slept through alarm

Child: Dry-night; child sticks on a sticker.

Parents: Tick where appropriate.

## Mattress covers

Mattress covers protect your bedlinen and keep your mattress dry if the child wets the bed. Further information: [www.rodger.nl](http://www.rodger.nl).

## Troubleshooting

a) System fails to give an alarm:

Replace the batteries, insuring that they are inserted in the right way.

Make sure that the contacts for the battery are clean.

Always use new batteries.

b) System continues to make noise:

The sensorclip is wet or dirty and requires cleaning.

## Help desk

For questions about the bedwetting treatment or the use of the bedwetting alarm, please get in contact with your supplier.

## Warranty

Rodger carries a six months warranty, after receipt of the device by the purchaser, against manufacturer's defects and against breakage or failure under normal use - at our discretion. This warranty does not cover batteries, battery lids, sensors, misuse or damage. We shall not be liable for more or other than the aforesaid warranty obligations.

This warranty is valid only if you are the first buyer and have the receipt. The warranty is not valid if the alarm is not used in accordance with the instructions or if repairs to the device are made by third parties. Rodger reserves the right to change specifications and warranty without prior notice.

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